

# ATISHAY JAIN

## EXPERIMENT 7: SODIUM BOROHYDRIDE REDUCTION OF BENZOIN

### PRE-LAB REPORT

DATE TO BE PERFORMED | 10/22/2019

TEACHING ASSISTANT: DR. PHIL SHIVOKEVICH

### PURPOSE :

The purpose of this experiment is to perform reduction of aldehyde/ ketone to form alcohol using metal hydride. For this experiment we will be using sodium borohydride ( $\text{NaBH}_4$ ) for reduction of BENZOIN.

### REAGENT TABLE:

COMPOUND	MW (g/mol )	AMOUN T NEEDED	MILLI- MOLES	BP	MP	DENSIT Y (g/ml)	HAZARDS
BENZOIN		0.5GM				-	
95% ETHANOL	-	4ML	-	-	-		
SODIUM BOROHYDRIDE		0.1GM	-			-	
ACETONE	-	JUST ENOUGH TO DISSOLVE	-		-	-	

HCl	-	0.3ml	-	-	-	-
ETHYL ACETATE	-	UPTO BASELINE (1CM)	-	-	-	-

## EXPERIMENTAL PROCEDURE

### Part I: Procedure for the Reduction of Benzoin

- Add benzoin (0.5 g) and ethanol (4 mL) to a 25-mL Erlenmeyer flask and swirl gently at room temperature for several minutes.
- Add sodium borohydride (0.1 g) using a micro-spatula in several small portions over 5 minutes.
- Gently swirl the mixture at room temperature for 20 additional minutes.
- Cool the mixture in an ice-water bath, then add water (5 mL) and 6M HCl (0.3 mL).
- Quenching: After 15 minutes, add an additional portion of water (2.5 mL).
- Collect the product via vacuum
- Allow the product to dry on the filter for 15 minutes, and then determine the crude yield and mp. Reserve a small portion (1-2 mg) of the crude material for TLC analysis.
- Recrystallize the crude solid using acetone.
- Determine the mp and yield (% and mass) of purified 1,2-diphenylethane-1,2-diol. Obtain an IR of the product with assistance from our TA.

### Part II: (TLC)

- You should dissolve a small amount of benzoin (the starting material), your recrystallized product and your reserved crude product (1-2 mg) in ethyl acetate in each of three vials.
- Spot 2 TLC plates, as shown below.
- Run the TLC plates in 9:1 CH<sub>2</sub>Cl<sub>2</sub>: ethanol
  - a. Add your eluant to a TLC developing chamber and carefully place the TLC plate in the developing chamber.
  - b. Remove the TLC plate, when solvent reaches about 1cm from the top, mark the solvent front and allow it to dry in your fume hood.
  - c. Once the plate is dry, visualize the spots under UV light and circle them with a pencil.